

Dual Frequency Ion Isolation in Quadrupole Ion Traps

Track Code: 2016-COOK-67503

Categories:

- Chemical Engineering
- Chemistry and Chemical Analysis

Keywords:

- Chemical Engineering
- Chemistry and Chemical Analysis
- Ion Trap
- Mass Spectrometry

Tandem mass spectrometry is a valuable tool in many fields, such as biological, chemical, and pharmaceutical research. The first step of this form of spectrometry is the isolation step when ions are separated by mass-to-charge ratio. The two most common isolation techniques are radio frequency/direct current (RF/DC) isolation and stored waveform inverse Fourier transform (SWIFT); however, both methods have drawbacks. In the former, the RF and DC must be electronically coupled together, adding unnecessary difficulty to the process. In the latter, excessively complex features are needed.

Researchers at Purdue University have developed a simpler ion isolation technique. This method is based on a single dipolar waveform with two frequency components. One frequency is chosen to eject ions lower in mass than the target ion, and the other frequency ejects ions higher in mass than the target ion. The dual frequency waveform isolates the targeted ion.

Advantages:

- Simple design
- Easy to implement

Potential Applications:

- Mass spectrometry
- Biological research
- Chemical research
- Pharmaceutical research
- Chemical analysis

People:

- Cooks, Robert Graham (Project leader)
- Snyder, Dalton Thomas

Intellectual Property:

Application Date: July 30, 2018

Type: NATL-Patent

Country of Filing: United States

Patent Number: 11,355,328

Issue Date: June 7, 2022

Application Date: April 13, 2017

Type: PCT-Patent

Country of Filing: WO

Patent Number: (None)

Issue Date: (None)

Application Date: April 13, 2016

Type: Provisional-Patent

Country of Filing: United States

Patent Number: (None)

Issue Date: (None)

Contact OTC:

Purdue Office of Technology Commercialization

The Convergence Center

101 Foundry Drive, Suite 2500

West Lafayette, IN 47906

Phone: (765) 588-3475

Fax: (765) 463-3486

Email: otcip@prf.org