

Continuous Lyophilization Roll-to-Roll Manufacturing and Metrology System

Track Code: 2017-CAKM-67672

Categories:

- Materials and Manufacturing

Keywords:

- Chemical Analysis
- Chemistry and Chemical Analysis
- Food and Nutrition
- Lyophilization
- Materials and Manufacturing
- Pharmaceuticals

Batch processing is a popular mechanism to freeze dry pharmaceuticals and food products. However, there is an increasing interest in the pharmaceutical industry to switch from batch to continuous freeze drying processes due to a need for pristine physical and chemical data regarding the process.

A Purdue University researcher has developed a new roll-to-roll freeze drying metrology system that addresses issues with current batch processes. This technology helps accelerate the process and quality control by optimizing the roll to roll freeze drying system. It can also access the physio-chemical changes that take place during freeze drying.

Advantages:

- Accelerates freeze drying process
- Increases data gathering

Potential Applications:

- Pharmaceuticals
- Food products
- Continuous freeze drying processes

People:

- Cakmak, Mukerrem (Project leader)

Intellectual Property:

Application Date: April 1, 2019

Type: NATL-Patent
Country of Filing: United States
Patent Number: 11,175,092
Issue Date: November 16, 2021

Application Date: November 15, 2021
Type: DIV-Patent
Country of Filing: United States
Patent Number: (None)
Issue Date: (None)

Application Date: October 10, 2017
Type: PCT-Patent
Country of Filing: WO
Patent Number: (None)
Issue Date: (None)

Application Date: October 10, 2016
Type: Provisional-Patent
Country of Filing: United States
Patent Number: (None)
Issue Date: (None)

Contact OTC:

Purdue Office of Technology Commercialization
The Convergence Center
101 Foundry Drive, Suite 2500
West Lafayette, IN 47906

Phone: (765) 588-3475
Fax: (765) 463-3486
Email: otcip@prf.org